

6/20 Sample C#2017

Start w/ 20 μ l DNA
26 μ l H₂O

coveris spin 4,200 rpm 1 min x 2

① DNA Repair

45 μ l DNA

8.5 μ l H₂O

6.5 FFPE repair buffer

2 μ l FFPE repair enzyme
62 μ l

in PCR tube and inc. at 20°C for 15 min

transfer to 1.5 ml tube

add 62 μ l of Ampure beads & inc RT 10 min
spin down & pellet on magnet
wash w/ 200 μ l 70% EtOH x 2

spin/remove last of EtOH & dry

resuspend w/ 48 μ l of H₂O inc 2 min

pellet & remove 46 μ l

Qubit 6.42 ng/ μ l (288.9 ng)

② End Prep

45 μ l DNA

7 μ l Ultra II end prep buff.

3 μ l " mix

5 μ l DCS

60 μ l

Continued on Page

Read and Understood By

Signed

Date

Signed

Date

6/20 Inc. 20°C 5min
USC 5min

transc to 1.5 ml tube
add 60µl AMPure beads

inc. 5min & flick

spindown & magnet

wash w/ 200µl 70% EtOH x2

spindown/remove rest & dry

resuspend w/ 31µl H₂O inc 2min

Qubit 8.64 ng/µl (259.2 ng)

③ Adapter Ligation

30µl DAA

20µl Adapter mix

50µl Blunt/TA Ligation mix
100µl

incubate 10 min.

④ Bead Binding

add 40µl Ampure Beads

flick & inc 5min

spindown & magnet

140µl ABB x2

15µl ELB

* Computer broken, refrigerated & run next day

Continued on

Read and Understood By

Signed

Date

Signed

Date